

# FACT SHEET

## The iBGStar® Blood Glucose Monitoring System



### About iBGStar®

iBGStar® is the first FDA cleared blood glucose meter that directly connects to the iPhone® and iPod touch®. With seamless integration into the lifestyles of people with diabetes, the real-time test results can be displayed, managed and downloaded on-the-go.

Results are presented on the Multi-Touch display of an iPhone® or iPod touch® via the iBGStar® Diabetes Manager App (free at the App Store<sup>SM</sup>), quickly after monitoring blood glucose. iBGStar® can also be used alone to measure blood glucose levels and results can be synchronized later to an iPhone® or iPod touch®.



### iBGStar® Technology

iBGStar® uses patented Dynamic Electrochemistry® technology that provides accurate results by helping to detect and correct for errors caused by differences in hematocrit levels, manufacturing variations and certain environmental conditions.

### iBGStar® Data Management

The iBGStar® Diabetes Manager App has a range of features that allow users multiple ways to analyze their glucose patterns over time. Consumers can record and track readings, carbohydrate intake and insulin doses (if taking insulin). Scorecards show individual test results and are color-coded so that high and low blood glucose results can be more easily identified. Trend charts with customizable views can be designed to determine if values are within a defined target range. The 'Share' function enables users to select specific data to send as an e-mail to their healthcare team.



In March 2010, Sanofi and AgaMatrix signed an agreement for the development, supply and commercialization of Blood Glucose Monitoring solutions. iBGStar® is the result of this agreement.

### Key iBGStar® Features

- Meets today's industry standards for accuracy
- Ultra compact size (width 0.94" x length 2.17" x height 0.39"); weight (9 g)
- Features a dot matrix organic light emitting diode (OLED) display that illuminates numbers, letters and symbols on a dark background, making them easily visible and eliminates the need for a backlight
- Storage for up to 300 test results with date and time on standalone iBGStar® meter; iPhone® and iPod touch® are able to store a nearly unlimited number of results depending on device capacity
- No coding required
- 24 hour customer support telephone hotline, 855-iBGStar (424-7827), and website support available upon commercial launch to patients and healthcare professionals

For iBGStar® product information and updates, please visit: [www.ibgstar.us](http://www.ibgstar.us)

Apple®, iPhone® and iPod touch® are trademarks of Apple Inc, registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Content purchased from the iTunes Store is for personal lawful use only. Don't steal music.

### Important Information

The iBGStar® meter and lancing device are for single patient individual use. Do not share them with anyone including other family members. Do not use on multiple patients. All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.



@Diabetes\_Sanofi



Facebook.com/sanofiUSDiabetes



discussdiabetes.com



Good Design 2011  
Medical



reddot design award

iBGStar® received the Good Design™ Award in 2011 for outstanding product design in the medical category from the Chicago Athenaeum of Architecture and Design and the European Centre for Architecture Art Design and Urban Studies.